



Recombinant Mouse Protein deltex-2 (Dtx2)

Product Code	CSB-MP854510MO
Abbreviation	Dtx2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R3P2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAMAPSSSLP QVYP SHVVVA VWEWQDGLGI WHPYSATVCS FIEQHFVRQR GQHFGGLGSLA HSIPLGQADP SLAPYIIDLP SWTQFRQNTG TMRSVRRHLF SQNSAPGQGI VWEWLGDDGS WVAYEARICD YLEQQVARGI QVVDLAPLGY NYTVNYATLT QTNKTSSFCR SVRRQVGPVY PVTSDIAVPR QMGLICFCQQ CLHGSGTG PV SGRYRHSMTN LPAYPAPQAP HRTTTVSGAH QAFAPYNKPS LSGARSAPRL NTTNPWAAAP PVAGNQSLFH SSLSHLGPQL LPSGPSTSSG ASAFPSPGPS SSSPGSAPTT VPVQMPKASR VQQALAGMTS VLSAIGLPVC LSRAPRPTGP PASRPASKSH SSVKRLRKMS VKEGAPKPEP EQVIRKYTEE LKVAPEEDCI ICMEKLAVAS GYS DMTDSKA LGPMVVGRLT KCSHAFHLLC LLAMYCNGNK DGSLQCPSCK TIYGEKTGTQ PWGKMEVFRF QMSLPGHEDC GTILIVYNIP HGIQGPEHPS PGKPFTARGF PRQCYLPDSP QGRKVLELLK VAWKRRIFT VGTSSTTGET DTVVWNEIHH KTEMDRNV TG HGYPDPNYLQ NVLAELAAQG VTEDCLEQQ
Source	Mammalian cell
Target Names	Dtx2
Protein Names	Recommended name: Protein deltex-2 Short name= Deltex2 Short name= mDTX2
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.