



# Recombinant Solanum tuberosum Alpha-1,4-glucan-protein synthase [UDP-forming] 2 (UPTG2)

<b>Product Code</b>	CSB-BP854630FIG
<b>Abbreviation</b>	UPTG2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8RU27
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Solanum tuberosum (Potato)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGSSVTPTP LLKDELDIVI PTIRNLDFLE MWRPFFQPYH LIIVQDGDPS KIINVPEGFD YELYNRNDIN RILGPKASCI SFKDSACRCF GYMVSKKKYI YTIDDDCFVA KDPSGKDINA LEQHIKLLC PSTPHFFNTL YDPYREGADF VRGYPFMRE GAATAVSHGL WLNIPDYDAP TQLVKPRERN TRYVDAVMTI PKGTLFPMCG MNLAFDRELI GPAMYFGLMG DGQPIGRYDD MWAGWCIKVI CDHLGLGVKT GLPYIWHKA SNPFVNLKKE YKGIYWQEEI IPFSQSATLP KDCTSVQQCY LELSKQVKEK LSTIDPYFTK LADAMVTWIE AWDELNPTGE GLAKLPSRTA PESRLH
<b>Source</b>	Baculovirus
<b>Target Names</b>	UPTG2
<b>Protein Names</b>	Recommended name: Alpha-1,4-glucan-protein synthase [UDP-forming] 2 EC=2.4.1.- Alternative name(s): Reversibly glycosylated polypeptide 2 Short name=RGP2 UDP-glucose:protein transglucosylase 2 Short name= UPTG 2
<b>Expression Region</b>	1-366
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.