



Recombinant *Oryza sativa* subsp. japonica Putative calcium-binding protein CML19 (CML19)

Product Code	CSB-EP854648OFG
Abbreviation	CML19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8RYJ8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVAATAEFRR VFSAFDRDAD GKISAAELRL CMKAALGEDM SAEAEALVS SADTDDDGLL DEEFETKLAV QLEMGDEEER CRGLMEAFRM YEMEGEGRIT PASLKRMLSK LGSHQGIEEC QTMICRFDLG GDGVISFEEF KIMMDA
Source	E.coli
Target Names	CML19
Protein Names	Recommended name: Putative calcium-binding protein CML19 Alternative name(s): Calmodulin-like protein 19
Expression Region	1-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.