



Recombinant Human GTPase RhebL1 (RHEBL1)

Product Code	CSB-MP855034HU
Abbreviation	RHEBL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TAI7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPLVRYRKVV ILGYRCVGKT SLAHQFVEGE FSEGYDPTVE NTYSKIVTLG KDEFHLHLVD TAGQDEYSIL PYSFIIGVHG YVLVYSVTSL HSFQVIESLY QKLHEGHGKT RVPVVLVGNK ADLSPEREVQ AVEGKLAES WGATFMESSA RENQLTQGIF TKVIQEIARV ENSYQGERRC
Source	Mammalian cell
Target Names	RHEBL1
Protein Names	Recommended name: GTPase RhebL1 Alternative name(s): Ras homolog enriched in brain like-1 c Short name= RhebL1c Ras homolog enriched in brain-like protein 1 Short name= Rheb-like protein 1 Rheb2
Expression Region	1-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.