



# Recombinant Human Centrosomal protein of 76 kDa (CEP76)

<b>Product Code</b>	CSB-YP855035HU
<b>Abbreviation</b>	CEP76
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TAP6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSLPPEKASE LKQLIHQQLS KMDVHGRIRE ILAETIREEL APDQQHLSTE DLIKALRRRG IIDDVMKELN FVTDSVEQEL PSSPKQPICF DRQSTLKKTN IDPTRRYLYL QVLGGKAFLE HLQEPEPLPG QVCSTFTLCL HYRNQRFRSK PVPCACEPDF HDGFLLEVHR ESLGDGTRMA DSTTMLSID PIHMVLIKTD IFGETTLVAS YFLEWRSVLG SENGVTSLTV ELMGVGTESK VSVGILNIKL EMYPLNQTL SQEVVNTQLA LERQKTAEKE RLFLVYAKQW WREYLQIRPS HNSRLVKIFA QDENGINRPV CSYVKPLRAG RLLDTPRQAA RFVNVLYGER APVIGGGGKQ EQWCTLLAFL CRNKGDCEDH ANLLCSLLLG YGLEAFVCVG TKAKGVPHAW VMTCGTDGAI TFWESLTGHR YIHKPTNPDE PPVAEQPKPL YPYRTIGCVF NHQMFLGNCQ PSDAVETCVF DLNDESKWKP MSEEAIKSVC APGATTSLPP FPPLCASTID ASVTSNEIEM QLRLLVSEHR KDLGLTTVWE DQLSYLLSPA LASYEFERTT SISAGNEEFQ DAIRRAVPDG HTFKGFPIHF VYRNARRAFA TCLRSPFCEE IICCRGDQVR LAVRVRVFTY PESACAVWIM FACKYRSVL
<b>Source</b>	Yeast
<b>Target Names</b>	CEP76
<b>Protein Names</b>	Recommended name: Centrosomal protein of 76 kDa Short name= Cep76
<b>Expression Region</b>	1-659
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.