



Recombinant Human SH3KBP1-binding protein 1 (SHKBP1)

| | |
|--------------------------|--|
| Product Code | CSB-MP855040HU |
| Abbreviation | SHKBP1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8TBC3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | <p>AAAATAAEG VPSRGPPGEV IHLNVGGKRF STSRQTLTWI PDSFFSLLS GRISTLKDET GAIFIDRDPT VFAPILNFLR TKELDPRGVH GSSLLHEAQF YGLTPLVRRL QLREELDRSS CGNVLFNGYL PPPVFPVKRR NRHSLVGPQQ LGGRPAPVRR SNTMPPNLGN AGLLGRMLDE KTPPSPSGQP EEPGMVRLVC GHHNWIAVAY TQFLVCYRLK EASGWQLVFS SPRLDWPIER LALTARVHGG ALGEHDKMVA AATGSEILLW ALQAEAGGSE IGVFHLGVPV EALFFVGNQL IATSHTGRIG VWNAVTKHWQ VQEVQPITSY DAAGSFLLLG CNNGSIYYVD VQKFPLRMKD NDLLVSELYR DPAEDGVTAL SVYLTPKTS D SGNWIEIAYG TSSGGVRVIV QHPETVGSGP QLFQFTFTVHR SPVTKIMLSE KHLISVCADN NHVRTWSVTR FRGMISTQPG STPLASFKIL ALESADGHGG CSAGNDIGPY GERDDQQVFI QKVVPSASQL FVRLSSTGQR VCSVRSVDGS PTTAFTVLEC EGSRRLLGSRP RRYLLTGQAN GSLAMWDLTT AMDGLGQAPA GGLTEQELME QLEHCELAPP APSAPSWGCL PPSPRISLT SLHSASSNTS LSGHRGSPSP PQAEARRRGG GSFVERCQEL VRSGPDLRRP PTPAPWPSSG LGTPLTPPKM KLNETSF</p> |
| Source | Mammalian cell |
| Target Names | SHKBP1 |
| Protein Names | Recommended name: SH3KBP1-binding protein 1 |
| Expression Region | 2-707 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final |



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.