



# Recombinant Human E3 ubiquitin-protein ligase DTX3L (DTX3L)

<b>Product Code</b>	CSB-BP855054HU
<b>Abbreviation</b>	DTX3L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TDB6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASHLRPPSP LLVRVYKSGP RVRKLESYF QSSKSSGGGE CTVSTQEHEA PGTFRVEFSE RAAKERVLLK GEHQILVDEK PVPIFLVPT NSIKKNTRPQ ISSLTQSQAE TPSGDMHQHE GHIPNAVDSC LQKIFLTVTA DLNCNLSKE QRAYITTLCPSIRKMEGHG IEKVCDFQD IERHQFLSE QFLESEQKQQ FSPSMTERKP LSQQERDSCI SPSEPETKAE QKSNYFEVPL PYFEYFKYIC PDKINSIEKR FGVNIEIQES SPNMVCLDFT SSRSGDLEAA RESFASEFQK NTEPLKQECV SLADSKQANK FKQELNHQFT KLLIKEKGGE LTLGTQDDI SAAKQKISEA FVKIPVKLFA ANYMMNVIEV DSAHYKLETT ELLQEISEIE KRYDICKSVS EKGQKTCILF ESKDRQVDLS VHAYASFIDA FQHASCQLMR EVLLKSLGK ERKHLHQTKF ADDFRKRHPN VHFVLNQESM TLTGLPNHLA KAKQYVLKGG GMSSLAGKKL KEGHETPMDI DSDDSKAASP PLKGSVSSEA SELDKKEKGI CVICMDTISN KKVLPCKKHE FCAPCINKAM SYKPICTCQ TSYGIQKGNQ PEGSMVFTVS RDSLPGYESF GTIVITYSMK AGIQTEEHPN PGKRYPGIQR TAYLPDNKEG RKVLKLLYRA FDQKLFTVG YSRVLGVSDV ITWNDIHHKT SRFGGPEMYG YPDPSTYKRV KEELKAKGIE
<b>Source</b>	Baculovirus
<b>Target Names</b>	DTX3L
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase DTX3L EC= 6.3.2.- Alternative name(s): B-lymphoma- and BAL-associated protein Protein deltex-3-like Rhysin-2 Short name= Rhysin2
<b>Expression Region</b>	2-740
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.