



Recombinant Human Glucosamine-6-phosphate isomerase 2 (GNPDA2)

Product Code	CSB-EP855058HU-B
Abbreviation	GNPDA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TDQ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRLVILDNYD LASEWAAKYI CNRIQFKPG QDRYFTLGLP TGSTPLGCIYK KLIEYHKNGH LSFKYVKTFN MDEYVGLPRN HPESYHSYMW NNFFKHIDID PNNAHILDGN AADLQAECDA FENKIKEAGG IDLFGGIGP DGHI AFNEPG SSLVSRTRLK TLAMDTILAN AKYFDGDLK VPTMALTGVG GTVMDAREVM ILITGAHKAF ALYKAIEEGV NHMWTVSAFQ QHPRTIFVCD EDATLELRVK TVKYFKGLMH VHNKLVDP LF SMK DGN
Source	E.coli
Target Names	GNPDA2
Protein Names	Recommended name: Glucosamine-6-phosphate isomerase 2 EC= 3.5.99.6 Alternative name(s): Glucosamine-6-phosphate deaminase 2 Short name= GNPDA 2 Short name= GlcN6P deaminase 2 Glucosamine-6-phosphate isomerase SB52
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.