



Recombinant Human Histone deacetylase complex subunit SAP25 (SAP25)

Product Code	CSB-MP855067HU
Abbreviation	SAP25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TEE9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTPLAPWDPK YEAKAGPRPV WGANCSGAS FSGRTLCHPS FWPLYEAASG RGLRPVAPAT GHWNGQQAPP DAGFPVVCCE DVFLSDPLL RGQRVPLYLS KAPQQMMGSL KLLPPPPIMS ARVLRPSPS RGPSTAWLSG PELIALTGLL QMSQGEPRPS SSAVGPPDHT SDPPSPCGSP SSSQGADLSL PQTPDTHCP
Source	Mammalian cell
Target Names	SAP25
Protein Names	Recommended name: Histone deacetylase complex subunit SAP25 Alternative name(s): 25 kDa Sin3-associated polypeptide Sin3 corepressor complex subunit SAP25
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.