



Recombinant Methanosarcina acetivorans Signal recognition 54 kDa protein (srp54)

Product Code	CSB-MP855097MFB
Abbreviation	srp54
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8THD0
Product Type	Recombinant Protein
Immunogen Species	Methanosarcina acetivorans (strain ATCC 35395 / DSM 2834 / JCM 12185 / C2A)
Purity	≥85% (SDS-PAGE)
Sequence	MVMEKLGDSL QGALKKLIGA GRIDERTVNE VVKDIQRALL QADVNVKLVLM GMSQRIKERA MKEEPAGMN PREHVIRIVY QELMEIIGKG AEIQLKPQTI MMVGLQGSGK TTSAAKLARY FQRKGLKAGV VAADTFRPGA YHQLKTLAEK LNVGFYGEEG NPDAVEITKN GLKALEKYDI RIVDTAGRHA LEADLIEEME RIHAVAKPDH KFMVLDAGIG QQASQQAHAH NDSVGITGVI ITKLDGTAKG GGALSAVSET KAPIAFIVGV ETPEDFEKFE ADRFISRLLG MGDCLKSLMEK AEESLSEEDV NVEALMQGRF TLKDMYKQLE AMNKMGPLKQ IMSMLPMGMG GMKFSDEMFO ATSDKMKNYK VIMDSMTEEE MTDPRVIGGS RIKRISKSGS CSSEVDRELL KYHKTMTAL KGFRGGKFNI QKMMKKKMMGM
Source	Mammalian cell
Target Names	srp54
Protein Names	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.