



Recombinant Arabidopsis thaliana Tubby-like F-box protein 3 (TULP3)

Product Code	CSB-EP855388DOA-B
Abbreviation	TULP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VY21
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSFKSLIQDM RGELGSISRK GFDVRFQYGR SRSQRVVQDT SVPVDAFKQS CWASMPPELL RDVLMRIEQS EDTWPSRKNV VSCAGVCRNW REIVKEIVRV PELSSKLTFP ISLKQPGPRG SLVQCYIMRN RSNQTYLYL GLNQAASNDD GKFLAAKRF RRPTCTDYII SLNCDDVSRG SNTYIGKLR NFLGTKFTVY DAQPTNPGTQ VTRTRSSRL SLKQVSPRIP SGNYPVAHIS YELNVLGSRG PRRMQCVMDA IPASAVEPGG TAPTQTELVH SNLDSFPSFS FFRSKSIRAE SLPSGPSSAA QKEGLLV LKN KAPRWHEQLQ CWCLNFNGRV TVASVKNFQL VAAPENGPAG PEHENVILQF GKVGKDVFTM DYQYPISAFQ AFTICLSSFD TKIACE
Source	E.coli
Target Names	TULP3
Protein Names	Recommended name: Tubby-like F-box protein 3 Short name= AtTLP3
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.