



Recombinant Arabidopsis thaliana COP9 signalosome complex subunit 6b (CSN6B)

Product Code	CSB-MP855408DOA
Abbreviation	CSN6B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W1P0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAPSSSSGLT FKLHPLVMLN ISDHFTRVKT QLNPPAASCA TGNGSNNADA MLLQNPRVYG CVIGLQRGRT VEIFNSFELI FDPALDTLDR SFLEKKQELY KKVFPDFYVL GWYSTGSDAT ESDMHIHKAL MDINESPVYV LLNPAINHAQ KDLPVTIYES EFHVIDGIPQ SIFVHTSYTI ETVEAERISV DHVAHLKPSD GGSAATQLAA HLTGIHSAIK MLNSRIRVLY QHIVAMQKGD KPCENSVLRQ VSSLLRSLPA AESEKFNENF LMEYNDKLLM SYLAMITNCT SNMNEVVDFK NTAYDKHSRR GGRTAFM
Source	Mammalian cell
Target Names	CSN6B
Protein Names	Recommended name: COP9 signalosome complex subunit 6b Short name= AtCSN6b Short name= Signalosome subunit 6b
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.