



Recombinant Arabidopsis thaliana Ubiquitin-like domain-containing CTD phosphatase (At4g06599)

Product Code	CSB-EP855411DOA
Abbreviation	At4g06599
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W3M6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASSSSSPTP SAAAAVSPLT EEELTLTVKW NGKEYTVRIC ADDSVAELKR RICLLTTVLP KRQKLLYPKI GNKLSDDSL LSSISFKPSL KMTMIGTVED DIIVDQAESP EIVDDFELGK EEAVDVKDKE VNKQKLRRRI DQYKINLRTP CRQGKLLVL DIDYTLFDHR STAENPLQLM RPYLHEFLTA AYA EYDIMIW SATSMKWVEL KMTGVLNN PNYKVTALLD HLAMITVQSD TRGIFDCKPL GLIWALLPEF YNPGNTIMFD DLRRNFVMNP QNGLTIKPFK KAHANRDTQD ELVKLTQYLL TIAELSDLSS LHHSRWESFS QDNVKKRRRQE
Source	E.coli
Target Names	At4g06599
Protein Names	Recommended name: Ubiquitin-like domain-containing CTD phosphatase EC=3.1.3.16 Alternative name(s): Nuclear proteasome inhibitor UBLCP1
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.