



Recombinant Arabidopsis thaliana COP9 signalosome complex subunit 3 (CSN3)

Product Code	CSB-EP855420DOA-B
Abbreviation	CSN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W575
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MIGAVNSVEA VITSIQGLSG SPEDLSALHD LLRGAQDSL R AEPGVNFSTL DQLDASKHSL GYLYFLEVLT CGPVSKEKAA YEIPIARFI NSCDAGQIRL ASYKFVSLCK ILKDHVIALG DPLRGVGPLL NAVQKLQVSS KRLTALHPDV LQLCLQAKSY KSGFSILSDD IVEIDQPRDF FLYSYYGGMI CIGLKRFQKA LELLYNVVT A PMHQVNAIAL EAYKKYILVS LIHNGQFTNT LPKCASTAAQ RSFKNYTGPY IELGNCYNDG KIGEALV V ARNAEFEEK NLGLVKQAVS SLYKRNILRL TQKYL TSLQ DIANMVQLGN AKEAEMHVLQ MIQDGQIHAL INQKDG MVR F LEDPEQYKSS EMIEIMDSVI QRTIGLSKNL LAMDESLSCD PLYLGKVGRE RQRYDFGDDF DTVPQKFSM
Source	E.coli
Target Names	CSN3
Protein Names	Recommended name: COP9 signalosome complex subunit 3 Short name= Signalosome subunit 3 Alternative name(s): Protein FUSCA 11
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.