



Recombinant Papiro hamadryas Tubulin beta-8 chain (TUBB8)

Product Code	CSB-BP855459ERJ
Abbreviation	TUBB8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WP13
Product Type	Recombinant Protein
Immunogen Species	Papiro hamadryas (Hamadryas baboon)
Purity	>85% (SDS-PAGE)
Sequence	MREIVLTQAG QCGNQIGAKF WEVISDEHAI DSAGTYHGDS RLQLERIDVY YNEACGGRYV PRAVLVDLEP GTLDSVRS GP FGQIFRPDSF IFGQCGAGNN WAKGHYTEGA ELMESVMDVV RKEAESCDCL QGFQLTHSLG GGTGSGMGTL LLSKIREEYP DRIINTFSIL PSPKVS DTVV EPYNATLSVH QLIENADETF CIDNEALYDI CSRTLKLPTP TYGDLNHLVS ATMSGVTTCL RFPGQLNADL RKLAVNMVPP PRLHFFMPGF APLTSRGSQQ YRALTVAE LT QQMFDARNMM AACDPRHGRY LTAAAI FRGR MPMREVDEQM FNIQKNSSY FADWLPHNVK TAVCDIPPRG LKMSATFIGN NTAIQELFKR VSEQFTAMFR RKAFLHWYTG EGMDEMEFTE AESNMNDLVS EYQQYQDATA EEEEFEEYAE EEEA
Source	Baculovirus
Target Names	TUBB8
Protein Names	Recommended name: Tubulin beta-8 chain
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.