



Recombinant Human Plasma cell-induced resident endoplasmic reticulum protein (PACAP)

Product Code	CSB-YP855472HU
Abbreviation	PACAP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WU39
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DRAPLTAT APQLDDEEMY SAHMPAHLRC DACRAVAYQM WQNLAKAETK LHTSNSGGRR ELSELVYTDV LDRSCSRNWQ DYGVREVDQV KRLTGPGGLSE GPEPSISVMV TGGPWPTLS RTCLHYLGEF GEDQIYEAHQ QGRGALEALL CGGPQGACSE KVSATREEL
Source	Yeast
Target Names	MZB1
Protein Names	Recommended name: Plasma cell-induced resident endoplasmic reticulum protein Short name= Plasma cell-induced resident ER protein Short name= pERp1 Alternative name(s): Proapoptotic caspase adapter protein
Expression Region	23-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.