



# Recombinant Human Protein spire homolog 2 (SPIRE2)

<b>Product Code</b>	CSB-MP855490HU
<b>Abbreviation</b>	SPIRE2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WWL2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MARAGSCGGA AAGAGRPEPW ELSLEEV LKA YEQPLNEEQ A WAVCFQGCRG LRGSPGRRLR DTGDLLL RGD GSVGAREPE A AEPATMVVPL ASSEAQTVQS LGFAIYRALD WGLDESEERE LSPQLERLID LMANNDSEDS GCGAADEGYG GPEEEEEAE G VPRSVRTFAQ AMRLCAARLT DPRGAQAHYQ AVCRALFVET LELRAFLARV REAKEMLQKL REDEPHLETP RAELDSLGH T DWARLWVQLM RELRRGVK LK KVQEQEFNPL PTEFQLTPFE MLMQDIRARN YKLRKVMVDG DIPPRVKKDA HELILDFIRS RPPLKQV SER RLRPLPPKQR SLHEKILEEI KQERRLRPVR GEGWAARGFG SLPCILNACS GDAKSTSCIN LSVTDAGGSA QRPRPRVLLK APTLAEMEEM NTSEEEESP C GEVTLKRDRS FSEHDLAQLR SEVASGLQSA THPPGGTEPP RPRAGSAHVW RPGSRDQGTC PASVSDPSHP LLSNRGSSGD RPEASMT PDA KHLWLEFSHP VESLALTVEE VMDVRRVLVK AEMEKFLQNK ELFSSLKKGK ICCCCRAKFP LFSWPPSCLF CKRAVCTSCS IKMKMPSKKF GHIPVYTLGF ESPQRVSAAK TAPIQRDIF QSLQGPQWQS VEEAFPHIYS HGCVLKDVCS ECTSFVADV V RSSRKSVDVL NTPRRSRQT QSLYIPNTRT LDFK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPIRE2
<b>Protein Names</b>	Recommended name: Protein spire homolog 2 Short name= Spir-2
<b>Expression Region</b>	1-714
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.