



Recombinant Human Ankyrin repeat and SOCS box protein 11 (ASB11)

Product Code	CSB-BP855501HU
Abbreviation	ASB11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WXH4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDGPFVYGF KNIFITMFAT FFFFKLLIKV FLALLTHFYI VKGNRKEAAR IAEEIYGGIS DCWADRSPLH EAAAQGRLLA LKTLIAQGVN VNLVTINRVS SLHEACLGGH VACAKALLEN GAHVNGVTVH GATPLFNACC SGSAACVNVL LEFGAKAQLV VHLASPIHEA VKRGHRECME ILLANNVNID HEVPQLGTPL YVACTYQRVD CVKKLLELGA SVDHGQWLDL PLHAAARQSN VEVHLLTDY GANLKRRNAQ GKSALDLAAP KSSVEQALLL REGPPALSQ LRLCVRKCLG RACHQAIHKL HLPEPLERFL LYQ
Source	Baculovirus
Target Names	ASB11
Protein Names	Recommended name: Ankyrin repeat and SOCS box protein 11 Short name= ASB-11
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.