



# Recombinant Human Dyslexia susceptibility 1 candidate gene 1 protein (DYX1C1)

<b>Product Code</b>	CSB-MP855504HU
<b>Abbreviation</b>	DYX1C1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WXU2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPLQVSDYSW QQTKTAVFLS LPLKGV CVRD TDVFCTENYL KVNFPFLFE AFlyAPIDDE SSKAKIGNDT IVFTLYKKEA AMWETLSVTG VDKEMMQRIR EKSILQAQER AKEATEAKAA AKREDQKYAL SVMMKIEEEE RKKIEDMKEN ERIKATKALE AWKEYQRKAE EQKKIQREEK LCQKEKQIKE ERKKIKYKSL TRNLASRNLA PKGRNSENIF TECLKEDSIP APRSVGSIKI NFTP RVFPTA LRESQVAEEE EWLHKQAEAR RAMNTDIAEL CDLKEEEKNP EWLKDKGNKL FATENYLA AI NAYNLAIRLN NKMP LLYLNR AACHLKLKLNL HKAIEDSSKA LELLMPPVTD NANARMKAHV RRGTAFCQLE LYVEGLQDYE AALKIDPSNK IVQIDAEKIR NVIQGTELKS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	DNAAF4
<b>Protein Names</b>	Recommended name: Dyslexia susceptibility 1 candidate gene 1 protein
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.