



Recombinant Human Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 2 (ARAP2), partial

Product Code	CSB-BP855511HU
Abbreviation	ARAP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WZ64
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	ARAP2
Protein Names	Recommended name: Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 2 Alternative name(s): Centaurin-delta-1 Short name=Cnt-d1 Protein PARX
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This protein contains ARF-GAP, RHO-GAP, ankyrin repeat, RAS-associating, and pleckstrin homology domains. The protein is a phosphatidylinositol (3,4,5)-trisphosphate-dependent Arf6 GAP that binds RhoA-GTP, but it lacks the predicted catalytic arginine in the RHO-GAP domain and does not have RHO-GAP activity. The protein associates with focal adhesions and functions downstream of RhoA to regulate focal adhesion dynamics.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.