



Recombinant Uncharacterized protein yghZ (yghZ)

Product Code	CSB-MP855528EOD
Abbreviation	yghZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X529
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MVWLANPERY GQMQYRYCGK SGLRLPALS LGLWHNFGHVN ALESQRAILR KAFDLGITHF DLANNYGPPP GSAEENFGRL LREDFAAAYRD ELIISTKAGY DMWPGPYGSG GSRKYLLASL DQSLKRMGLE YVDIFYSHRV DENTPMEETA SALAHAVQSG KALYVGISSY SPERTQKMVE LLHEWKIPLL IHQPSYNLLN RWVDKSGLLD TLQNNGVGC I AFTPLAQGLL TGKYLNGIPE DSRMHREGNK VRGLTPKMLT EANLNSLRLL NEMAQQRGQS MAQMALSPLL KDERVTSVLV GASRAEQLEE NVQALNNLTF STEELAQIDQ HIADGELNLW QASSDK
Source	Mammalian cell
Target Names	gpr
Protein Names	Recommended name: Uncharacterized protein yghZ
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.