



Recombinant 2,3,4,5-tetrahydropyridine-2,6-dicarboxylate N-succinyltransferase (dapD)

Product Code	CSB-MP855559EOD
Abbreviation	dapD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X8Y7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MQQLQNIET AFERRAEITP ANADTVTREA VNQVIALLDLS GALRVAEKID GQWVTHQWLK KAVLLSFRIN DNQVIEGAES RYFDKVPKMF ADYDEARFQK EGFRVVPAA VRQGAFIARN TVLMPSYVNI GAYVDEGTMV DTWATVGSCA QIGKNVHLSG GVGIGGVLEP LQANPTIIED NCFIGARSEV VEGVIVEEGS VISMGVYIGQ STRIYDRETG DIHYGRVPAG SVVVSIGNLPS KDGKYSLYCA VIVKKVDAKT RGKVGINELL RTID
Source	Mammalian cell
Target Names	dapD
Protein Names	Recommended name: 2,3,4,5-tetrahydropyridine-2,6-dicarboxylate N-succinyltransferase EC= 2.3.1.117 Alternative name(s): Tetrahydrodipicolinate N-succinyltransferase Short name= THDP succinyltransferase Short name= THP succinyl
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.