



Recombinant 1-deoxy-D-xylulose-5-phosphate synthase (dxs)

Product Code	CSB-EP855606EOD
Abbreviation	dxs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8XE76
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	SFDIAKYPT LALVDSTQEL RLLPKESLPK LCDELRRYLL DSVSRSSGHF ASGLGTVELT VALHYVYNTF FDQLIWDVGH QAYPHKILTG RRDKIGTIRQ KGGLHPFPWR GESEYDVLSV GHSSTSISAG IGIAVAAEKE GKNRRTVCVI GDGAITAGMA FEAMNHAGDI RPDMLVVLND NEMSIENVG ALNNHLAQLL SGKLYSSLRE GGKKVFSGVP PIKELLKRTE EHIKGMVVPGLFEELGFNY IGPVDGHDVL GLITTLKNMR DLKGPQFLHI MTKKGRGYEP AEKDPITFHA VPKFDPSSGC LPKSSGGLPS YSKIFGDWLC ETAAKDNKLM AITPAMREGS GMVEFSRKFP DRYFDVAIAE QHAVTFAAGL AIGGYKPIVA IYSTFLQRAY DQVLHDVAIQ KLPVLFVIDR AGIVGADGQT HQGAFDLSYL RCIPEMVIMT PSDENECRQM LYTGYYHNDG PSAVRYPRGN AVGVELTPLE KLPIGKGIVK RRGEKLAILN FGTLMPEAAK VAESLNATLV DMRFVKPLDE TLILEMAASH EALVTVEENA IMGGAGSGVN EVLMAHRKPV PVLNIGLPDF FIPQGTQEEM RAELGLDAAG MEAKIKAWLA
Source	E.coli
Target Names	dxs
Protein Names	Recommended name: 1-deoxy-D-xylulose-5-phosphate synthase EC= 2.2.1.7 Alternative name(s): 1-deoxyxylulose-5-phosphate synthase Short name= DXP synthase Short name= DXPS
Expression Region	2-620
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.