



Recombinant Tetraodon fluviatilis Proto-oncogene c-Fos (fos)

Product Code	CSB-EP856231TLV-B
Abbreviation	fos
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91496
Product Type	Recombinant Protein
Immunogen Species	Tetraodon fluviatilis (Green pufferfish) (Chelonodon fluviatilis)
Purity	>85% (SDS-PAGE)
Sequence	MMFTSFNAEC DSSSRCSASP SDNVYYPSPA GSYSSMGSPQ SQDLTDLTAS SASFVPTVTA ISTSPDLQWM VQPLVSSVAP SRRHPYSPS PSYKRTVMRS GASKPHAKRG RVEQTTPEEE EKKRIRRERN KQAAAKCRNR RRELTDSLQA ETDQLEAEKS SLQNDIANLL KEKERLEFIL AAHQPIKIP SQMDSDFPVV SMSPVHAYLS TAASTQPQTS VPEATTVTSS HSTFTSTSNS IFGSNSDSL STATVSDSVV KMTDLESSVL EESLDLLAKT EVETVEVPDV NLSSSLYTAQ DWEPLHATIG SSDFEPLCTP VVTCTPACTT ITSSFVFTFP EAETFPTCCV AHRGNSNSND QSSDSLSSPT LLAL
Source	E.coli
Target Names	fos
Protein Names	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.