



Recombinant Maize streak virus genotype E Capsid protein (V1)

Product Code	CSB-EP856285MIAM
Abbreviation	V1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91MG3
Product Type	Recombinant Protein
Immunogen Species	Maize streak virus genotype E (isolate Pat) (MSV)
Purity	>85% (SDS-PAGE)
Sequence	MSTSKRKRAD EVQWNKRSTK KKASAPPVKK TGGKADRPSL QIQTLLHSGD TMITVPSGGG CDLISTYARG SDEGNRHTSE TLTYKVGVDY HFVADAGACR YSNRGTGVMW LVYDTPGGN SPSTKDIFAY PDALVAWPTT WKVSRELCHR FVVKRRWLFT METDGRIGSD IPPTNQSWPP CKRNIDFHKF TSGLGVRTQW KNVTDGGVGA IQRGALYMVI APGNVTFFTA HGQTRLYFKS VGNQ
Source	E.coli
Target Names	V1
Protein Names	Recommended name: Capsid protein Alternative name(s): Coat protein Short name= CP
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.