



Recombinant Mouse LIM domain-containing protein ajuba (Ajuba)

Product Code	CSB-MP856315MO
Abbreviation	Ajuba
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91XC0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MERLGEKASR LLEKLRLSDS GSAKFGRRKG EASRSGSDGT PGAGKGRLSG LGGPRKSGHR GANGGPGDEP LEPAREQGGL DAERNARGSF EAQRFEGSFP GGPPPTRALP LPLSSPPDFR LETTAPALSP RSSFASSAS DASKPSSPRG SLLLDGAGAS GAGGSRPCSN RTSGISMGYD QRHGSPLPAG PCLFGLPLTT APAGYPPGAP SAYPELHAAL DRLCAHRSVG FGCQESRHSY PPALGSPGAL TGAVVGTAGP LERRGAQPGR HSVTGYGDCA AGARYQDELTA ALLRLTVATG GREAGARGEP SGIEPSGLEE SPGPFVPEAS RSRIREPEAR EDYFGTCIKC NKGIYQSNAC QALDSLYHT QCFVCCSCGR TLRCKAFYSV NGSVYCEEDY LFSGFQEAAE KCCVCGHLIL EKILQAMGKS YHPGCFRCIV CNKCLDGVPF TVDFSNQVYC VTDYHKNYAP KCAACGQPIL PSEGCEDIVR VISMDRDIYHF ECIYHCEDCRM QLSDEEGCCC FPLDGHLLCH GCHMQRLSAR QPSTNYI
Source	Mammalian cell
Target Names	Ajuba
Protein Names	Recommended name: LIM domain-containing protein ajuba
Expression Region	1-547
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.