



Recombinant Mouse Sentrin-specific protease 2 (Senp2)

Product Code	CSB-EP856342MO-B
Abbreviation	Senp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91ZX6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MYRWLAKVLG TILRLCERPA PGARALLKRR RSSSTLFSTA VDTDEIPAKR PRLDCFIHQV KNSLYNAASL FGFPFQLTTK PMVSSACNGT RNVAPSGEVF SNSSSCELMS SGSCSSMLKL GNKSPNGISD YPKIRVTVTR DQPRRVLPSF GFTLKSEGYN RRPSGRRHSK SNPESLTKW PQEQGVTEMI SEEGGKGVRR PHCTVEEGVQ KDEREKYRKL LERLKEGAHG STFPPTVSHH SSQRIQMDTL KTKGWVEEQN HGVRTTHFVP KQYRVVETRG PLCSMRSEKR YSKGKADTEK VVGLRFEKEG TRGHQMEPDL SEEVSARLRL GSGSNGLLRR KISVLEIKEK NFPSKEKDRR TEDLFEFTED MEKEISNALG HGPPDEILSS AFKLRITRGD IQTLKNYHWL NDEVINFYMN LLVERSKKQG YPALHAFSTF FYPKLGSGGY QAVKRWTKGV NLFEQELVLV PIHRKVHWSL VVMDLRKKCL KYLDMSGQKG HRICEILLQY LQDESKTKRN TDLNLEWTH YSMKPHEIPQ QLNGSDCGMF TCKYADYISR DKPITFTQHQ MPLFRKKMVW EILHQQLL
Source	E.coli
Target Names	Senp2
Protein Names	Recommended name: Sentrin-specific protease 2 EC= 3.4.22.68 Alternative name(s): Axa2 SUMO-1 protease 1 Short name= SuPr-1 SUMO-1/Smt3-specific isopeptidase 2 Short name= Smt3ip2 Sentrin/SUMO-specific protease SEN
Expression Region	1-588
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.