



Recombinant *Xenopus laevis* Sonic hedgehog protein (shh)

Product Code	CSB-MP856343XBE
Abbreviation	shh
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92000
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	CGPGRG IGKRRHPKKL TPLAYKQFIP NVAEKTGAS GRYEGKITRN SDCFKELTPN YNPDIMFKDE ESTGADRLMT QRCKDKLNAL AISVMNQWPG VKLRVTEGWD EDGHHLEESL HYEGRAVDIT TSDRDRSKYG MLGRLAVEAG FDWVYYESKA HIHCSVKAEN SVAAKSGGCF PAGARVMVEF GGTKAVKDLR PGDRVLSDDP QGNLLYSDFL MFIDQERDVK KLFYVIETSQ RKIRLTA AHL LFVAQTKVNG TRSFKSVFAS NIQPGDLIYT ADPKTMTLKA VKVEKVDLEE DTGAYAPLTA HGTVIDQVL ASCYAVIEEH TWAHLAFAPL RFGMSLSSYI YPRDSSPPSG LQPHHQVDLQ SHHQVDLQSH HQVDLQSHHQ LEGIHWYSQL LYQIGTWLLD SNSLHPLGMA TKSS
Source	Mammalian cell
Target Names	shh
Protein Names	Recommended name: Sonic hedgehog protein Alternative name(s): VHH-1 X-SHH Cleaved into the following 2 chains: 1. Sonic hedgehog protein N-product 2. Sonic hedgehog protein C-product
Expression Region	25-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.