



Recombinant Mouse Serotransferrin (Tf)

Product Code	CSB-MP856363MO
Abbreviation	Tf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92111
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	V PDKTVKWCAV SEHENTKCIS FRDHMKTVLP PDGPRRLACVK KTSYPDCIKA ISASEADAMT LDGGWVYDAG LTPNNLKPVA AEFYGSVEHP QTYYYAVAVV KKGTFDFQLNQ LEGKKSCHTG LGRSAGWVIP IGLLFCKLSE PRSPLEKAVS SFFSGSCVPC ADPVAFPKLC QLCPGCGCSS TQPFFGYVGA FKCLKDGGGD VAFVKHTTIF EVLPEKADRQ QYELLCLDNT RKPVDQYEDC YLARIPSHAV VARKNNGKED LIWEILKVAQ EHFVGKSKD FQLFSSPLGK DLLFKDSAFG LLRVPPRMDY RLYLGHNYVT AIRNQQEGVC PEGSIDNSPV KWCALSHLER TKCDEWSIIS EGKIECESAE TTEDCIEKIV NGEADAMTLD GGHAYIAGQC GLVPVMAEYY ESSNCAIPSQ QGIFPKGYA VAVVKASDTS ITWNNLKGKK SCHTGVDRTA GWNIPMGMLY NRINHCKFDE FFSQGCAPGY EKNSTLCDLC IGPLKCAPNN KEEYNGYTGA FRCLVEKGDV AFVKHQTVLD NTEGKNPAEW AKNLKQEDFE LLCPDGTRKP VKDFASCHLA QAPNHVVVSR KEKAARVKAV LTSQETLFGG SDCTGNFCLF KSTTKDLLFR DDTKCFVKLP EGTTPKEYLG AEYMQSVGNM RKCSTSRLL EACTFHKH
Source	Mammalian cell
Target Names	Tf
Protein Names	Recommended name: Serotransferrin Short name= Transferrin Alternative name(s): Beta-1 metal-binding globulin Siderophilin
Expression Region	20-697
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a glycoprotein with an approximate molecular weight of 76.5 kDa. It is thought to have been created as a result of an ancient gene duplication event that led to generation of homologous C and N-terminal domains each of which binds one ion of ferric iron. The function of this protein is to transport iron from the intestine, reticuloendothelial system, and liver parenchymal cells to all proliferating cells in the body. This protein may also



have a physiologic role as granulocyte/pollen-binding protein (GPBP) involved in the removal of certain organic matter and allergens from serum.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.