



Recombinant Human Protein FAM3C (FAM3C)

Product Code	CSB-MP856401HU
Abbreviation	FAM3C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92520
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QVFEIK MDASLGNLFA RSALDTAARS TKPPRYKCGI SKACPEKHFA FKMASGAANV VGPKICLEDN VLMSGVKNNV GRGINVALAN GKTGEVLDTK YFDMWGGDVA PFIEFLKAIQ DGTIVLMGTY DDGATKLNDE ARRLIADLGS TSITNLGFRD NWFVFCGGKGI KTKSPFEQHI KNNKDTNKYE GWPEVVEMEG CIPQKQD
Source	Mammalian cell
Target Names	FAM3C
Protein Names	Recommended name: Protein FAM3C Alternative name(s): Interleukin-like EMT inducer
Expression Region	25-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene is a member of the family with sequence similarity 3 (FAM3) family and encodes a secreted protein with a GG domain. A change in expression of this protein has been noted in pancreatic cancer-derived cells. Alternate transcriptional splice variants which encode the same protein have been characterized.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.