



Recombinant Rat SH3 domain-containing kinase-binding protein 1 (Sh3kbp1)

Product Code	CSB-MP856412RA
Abbreviation	Sh3kbp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q925Q9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MVEAIVEFDY QAQHDELTI SVGEVITNIR KEDGGWWEQG INGRRLFPD NFKVREIKKDV KKDLLSNKAP EKPMHDVSSG NSLLSSETIL RTNKRGERRR RRCQVAFSYL PQNDDELELK VGDIIEVVGE VEEGWWEQVL NGKTGMFPSN FIKELSGESD ELGISQDEQL SKSRPEGFLP ASLLPFAHG AKGKTTFEGT ILYRAAPGKT EGHRRYYSR ETTGSESDGG DSSSTKSEGA NGTVATAAIQ PKKVKGVGFG DIFKDKPIKL RPRSIEVEND FLPVEKTIGK KLPPATSTPD PSKTEMDSRT KTKDYCKVIF PYEAQNDEL TIKEGDIVTL INKDCIDVGW WEGELNGRRG VFPDNFVKLL PSDFDKEGGR PKKPPPSAP VIKQGAGTTE RKHEIKKIPP ERPETLPNRT EEKERPEREP KLDLQKPSVP AIPPKKPRPP KTNSLNRPGV LPPRRPERPV GPLTHTRGDS SKIDLAGSTL SGILDKDLSD RSNDIDLEGF DSVISSTEKL SHPTTSRPAK TGRRPPSQSL TSSSLSSPDI FDSPSPEEDK EEHISLAHRG IDVSKKTSRT VTISQVSDNK ASLPPKPGTM AAASSGPASL SSVASSPMSS SLGTAGQRAS SPSLFAEGK AKTESAVSSQ AAIEELKMQV RELRTIETM KDQKREIKQ LLSLDEEKK IRLRLQMEVN DIKKALQSK
Source	Mammalian cell
Target Names	Sh3kbp1
Protein Names	Recommended name: SH3 domain-containing kinase-binding protein 1 Alternative name(s): Regulator of ubiquitous kinase Short name= Ruk SH3-containing, expressed in tumorigenic astrocytes
Expression Region	1-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.