



Recombinant Human Hermansky-Pudlak syndrome 1 protein (HPS1)

Product Code	CSB-MP856438HU
Abbreviation	HPS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92902
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKCVLVATEGAEVLFYWTDQEFEEESLRLKFGQSENEEEELPALEDQLSTLLAP VISSMTMLEKLSDTYTCFSTENGNFLYVLHLFGCELFIAINGDHTESEGDRLRK LYVLKYLFEVHFGLVTVDGHLIRKELRPPDLAQRVQLWEHFQSLLWYTSRLRE QEQCFAVEALERLIHPQLCELCIEALERHVIQAVNTSPERGGEEALHAFLLVHS KLLAFYSSHSASSLRPADLLALILLVQDLYPSESTAEDDIQPSPRRARSSQNIPV QQAWSPHSTGPTGGSSAETETDSFSLPEEYFTPAPSPGDQSSGEDRRKAGG NNS
Source	Mammalian cell
Target Names	HPS1
Protein Names	Recommended name: Hermansky-Pudlak syndrome 1 protein
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length of Isoform III
Target Details	This gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. The encoded protein is a component of three different protein complexes termed biogenesis of lysosome-related organelles complex (BLOC)-3, BLOC4, and BLOC5. Mutations in this gene are associated with Hermansky-Pudlak syndrome type 1. Multiple transcript variants encoding distinct isoforms have been identified for this gene; the full-length sequences of some of these have not been determined yet.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.