



Recombinant Human Lipoma-preferred partner (LPP)

Product Code	CSB-BP856611HU
Abbreviation	LPP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93052
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSHPSWLPPK STGEPLGHVP ARMETTHSFG NPSISVSTQQ PPKKFAPVVA PKPKYNPYKQ PGGEGDFLPP PPPPLDDSSA LPSISGNFPP PPPLDEEAFK VQGNPGGKTL EERRSSLDAE IDSLTSILAD LECSSPYKPR PPQSSTGSTA SPPVSTPVTG HKRMVIPNQP PLTATKKSTL KPQPAPQAGP IPVAPIGTLK PQPQPVPASY TTASTSSRPT FNVQVKSAP SPHYMAAPSS GQIYGSGPQG YNTQPVPVSG QCPPPSTRGG MDYAYIPPPG LQPEPGYGYA PNQGRYYEGY YAAGPGYGGR NDSPTYGQQ GHPNTWKREP GYTPPGAGNQ NPPGMYPVTG PKKTYITDPV SPCAPPLQP KGGHSGQLGP SSVAPSRPE DELEHLTKKM LYDMENPPAD EYFGRCARCG ENVVGEGETGC TAMDQVFHVD CFTCIICNNK LRGQPFYAVE KKAYCEPCYI NTLEQCNVCS KPIMERILRA TGKAYHPHCF TCVMCHRSLD GIPFTVDAGG LIHCIEDFHK KFAPRCSVCK EPIMPAPGQE ETVRIVALDR DFHVHCYRCE DCGGLLSEGD NQGCYPLDGH ILCKTCNSAR IRVLTAKAST DL
Source	Baculovirus
Target Names	LPP
Protein Names	Recommended name: Lipoma-preferred partner Alternative name(s): LIM domain-containing preferred translocation partner in lipoma
Expression Region	1-612
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.