



# Recombinant Human Lipoma-preferred partner (LPP)

<b>Product Code</b>	CSB-EP856611HU
<b>Abbreviation</b>	LPP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q93052
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSHPSWLPPK STGEPLGHVP ARMETTHSFG NPSISVSTQQ PPKKFAPVVA PKPKYNPYKQ PGGEGDFLPP PPPPLDDSSA LPSISGNFPP PPPLDEEAFK VQGNPGGKTL EERRSSLDAE IDSLTSILAD LECSSPYKPR PPQSSTGSTA SPPVSTPVTG HKRMVIPNQP PLTATKKSTL KPQPAPQAGP IPVAPIGTLK PQPQPVPASY TTASTSSRPT FNVQVKSAP SPHYMAAPSS GQIYGSGPQG YNTQVPVSG QCPPPSTRGG MDYAYIPPPG LQPEPGYGYA PNQGRYYEGY YAAGPGYGGR NDSPTYGQQ GHPNTWKREP GYTPPGAGNQ NPPGMYPVTG PKKTYITDPV SAPCAPPLQP KGGHSGQLGP SSVAPSRPE DELEHLTKKM LYDMENPPAD EYFGRCARCG ENVVGEGETGC TAMDQVFHVD CFTCIICNNK LRGQPFYAVE KKAYCEPCYI NTLEQCNVCS KPIMERILRA TGKAYHPHCF TCVMCHRSLD GIPFTVDAGG LIHCIEDFHK KFAPRCVCK EPIMPAPGQE ETVRIVALDR DFHVHCYRCE DCGGLLSEGD NQGCYPLDGH ILCKTCNSAR IRVLTAKAST DL
<b>Source</b>	E.coli
<b>Target Names</b>	LPP
<b>Protein Names</b>	Recommended name: Lipoma-preferred partner Alternative name(s): LIM domain-containing preferred translocation partner in lipoma
<b>Expression Region</b>	1-612
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.