



Recombinant Uncharacterized oxidoreductase F53C11.3 (F53C11.3)

Product Code	CSB-EP856631CXY-B
Abbreviation	F53C11.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93761
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAHCLHPEKF FPVVKSVLALP PGSLNGKVAL VTGGGTGLGK AIATTF AHLG ASVAIAARRL DVLEKTADEI RSSTGGVCEP FQMDVKDPAK VAKAFDAVEK KLGHTPDILI NNAAGNFIMA TERLSPNAYG TIIDIVLKGT LHVTTTELGRRCIQQKRGASV LSITTLYAQS GAPFVVPSAV SKAGVENMTK SLASEWAKHG LRFNAIAPGP IPTEGAFGRL FAGELKDSGD AMKASVPVGR LGHPPEIANL AAFMSSDFMS WMNGAIDFD GGQQHIHHGS HMGQFLHEWD NEKWEETENL IRGRTGKEKK SKL
Source	E.coli
Target Names	F53C11.3
Protein Names	Recommended name: Uncharacterized oxidoreductase F53C11.3 EC= 1.-.-.-
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.