



Recombinant Arabidopsis thaliana 40S ribosomal protein S4-1 (RPS4A)

Product Code	CSB-MP856664DOA
Abbreviation	RPS4A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93VH9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MARGLKKHLK RLNAPKHWML DKLGGAFAPK PSSGPHKSRE CLPLVLIIRN RLKYALTYRE VISILMQRHI QVDGKVRTDK TYPAGFMDVV SIPKTNENFR LLYDTKGRFR LHSIKDEEAK FKLCKVRSIQ FGQKGIPYLN TYDGRTIRYP DPLIKPNDTI KLDLEENKIV EFIKFDVGNV VMVTGGRNRG RVGVIKNREK HKGSFETIHI QDSTGHEFAT RLGNYVTIGK GTKPWVSLPK GKGIKLTIIE EARKRLSAQQ A
Source	Mammalian cell
Target Names	RPS4A
Protein Names	Recommended name: 40S ribosomal protein S4-1
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.