



Recombinant Cell death specification protein 2 (ces-2)

Product Code	CSB-MP856690CXY
Abbreviation	ces-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94126
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MDFHRALSAL FTNQA AVQPL LGSLGF PFND GTSILTTALA AQSGGKKLDT PLGILPFDSL PTTNLLTPTK KIKLEDELCA SPVSSRSSTV SSSH FSSPQR SPSRKMSVPI PEEKKDSAYF ERRRKNNDAA KR SRDARRQK EEQIASKAHA LERENMQLRG KVSSLEQEAA QLRFL FFSKI SPANSEVSCE SNDSTETNDS NDSKSDSTIE V
Source	Mammalian cell
Target Names	ces-2
Protein Names	Recommended name: Cell death specification protein 2
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.