



Recombinant Protein maternal effect lethal 26 (mel-26)

Product Code	CSB-BP856698CXY
Abbreviation	mel-26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94420
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MEPRIDGGVF IGGIGNSGNE MCSNGVPALG VSSQTEIKVE KVQHTWTVKN FSHCYQEYLE NFVYLQRGDE QLTWSIKIYP KNGENNKDF VFLCLNRVIN NNVKAGKIGF KSQFKLRTAE NKDIEMRIHP NPSHSDYVSY IKRDVLF PQI MPRDMIIVNV EIDVAVETIT TTNEPIQFEP TNSEQQLIED YQRLFSQELL CDFAINVNGK IIRAHKAVLA ARSPVFNAML THQDTDEAKS SMMYINDMDY DVIYEMVYYI YCGRCNKDIT DMATALLIAA DKYRLEELKS HCEKYLVENI NIENACSLLI IGDLYSAPKL RKRAVTYILA RPKNVTGTPG WEDILKGHPN LITDIFSQID RQSSTGATSS VSNLPGVPM D IPGITGNIVP PPSGL
Source	Baculovirus
Target Names	mel-26
Protein Names	Recommended name: Protein maternal effect lethal 26
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.