



Recombinant Arabidopsis thaliana F-box protein PP2-B11 (PP2B11)

Product Code	CSB-EP856716DOA-B
Abbreviation	PP2B11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q949S5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNNLPEDCIA KILSLTTPLD VCRLSAVSSI FRSAAGSDDV WNHFLPADFP AGFAAPAGLP TRKQLFFSLV DNPLLINGTL LSFSLERKSG NKCYYMMAARA LNIVWGHEQR YWHWISLPNT RFGEVAELIM VWWLEITGKI NITLLSDDTL YAAAYFVFKWN HSPYGFQRQPV ETSLVLADTE STDNVVQPSM ISLMQDSGGE EGQSPVLRD GWYEVELGQF FKRRGDLGEI EMSLKETKGP YEKKGLIVYG IEIRPVP
Source	E.coli
Target Names	PP2B11
Protein Names	Recommended name: F-box protein PP2-B11 Alternative name(s): Protein PHLOEM PROTEIN 2-LIKE B11 Short name= AtPP2-B11 SKP1-interacting partner 12
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.