



Recombinant Arabidopsis thaliana BTB/POZ and TAZ domain-containing protein 2 (BT2)

Product Code	CSB-BP856726DOA
Abbreviation	BT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94BN0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEAVLVAMSV PATTEDDGFS LITDKLSYNL TPTSDVEIVT SDNRRIPAHS GVLASASPVL MNIMKKPMRR YRGCGSKRVI KILGVPCDAV SVFIKFLYSS SLTEDEMERY GIHLLALSHV YMVTQLKQRC SKGVVQRLTT ENVVDVLQLA RLCDAPDVCL RSMRLIHSQF KTVEQTEGWK FIQEHDPFLE LDILQFIDDA ESRKKRRRRH RKEQDLYMQL SEAMECIEHI CTQGCTLVGP SNVVDNNKKS MTAEKSEPCK AFSTCYGLQL LIRHFAVCKR RNNDKGCLRC KRMLQLFRLH SLICDQPDSC RVPLCRQFRK RGEQDKKMGE DTKWKLLVTR VVSAKAMTSL CQSKKNKCEQ AQGV
Source	Baculovirus
Target Names	BT2
Protein Names	Recommended name: BTB/POZ and TAZ domain-containing protein 2 Alternative name(s): BTB and TAZ domain protein 2
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.