



# Recombinant Arabidopsis thaliana NAC domain-containing protein 89 (NAC089)

<b>Product Code</b>	CSB-EP856735DOA
<b>Abbreviation</b>	NAC089
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q94F58
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDTKAVGVSK DTAASMEAST VFPGFKFSPT DVELISYYLK RKMDGLERSV EVIPDLEIYN FEPWDLDPDKS IVKSDSEWFF FCARGKKYPH GSQNRRTKM GYWKATGKER DVKSGSEVIG TKRTL VFHIG RAPKERTDW IMHEYCVKGV SLDDAMVVCR VRRNKEYNSG TSQKAPKPN SAEKHAKVQN GATSSGSPSD WDNLVDFYLA GESGEKLLAE MAESSEN LQV DNDEFFADI LRDEIINLDE AVMTGNT PNE VPTLESASME IRVLP LNMI DKQMSS LLEE RPSQKKKGKD ATESSLSCFV GLYSIKSVNK ARWDVIIGVV ALIAMLFYLE
<b>Source</b>	E.coli
<b>Target Names</b>	NAC089
<b>Protein Names</b>	Recommended name: NAC domain-containing protein 89 Short name= ANAC089 Alternative name(s): Protein FRUCTOSE-SENSING QUANTITATIVE TRAIT LOCUS 6
<b>Expression Region</b>	1-340
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.