



Recombinant Bubalus bubalis Interleukin-2 (IL2)

Product Code	CSB-BP856812BWX
Abbreviation	IL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95KP3
Product Type	Recombinant Protein
Immunogen Species	Bubalus bubalis (Domestic water buffalo)
Purity	>85% (SDS-PAGE)
Sequence	APTSSSTGNT MKEVKSLLLD LQLLLEKVKN PENLKLSRMH TFNFYVPKVN ATELKHLKCL LEELKLLLEV LNLAPSKNLN PREIKDSMDN IKRIVLELQG SETGFTCEYD DATVKA VEFL NKWITFCQSI YSTMT
Source	Baculovirus
Target Names	IL2
Protein Names	Recommended name: Interleukin-2 Short name= IL-2 Alternative name(s): T-cell growth factor Short name= TCGF
Expression Region	21-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.