



Recombinant *Drosophila melanogaster* N-alpha-acetyltransferase 30 (CG11412)

Product Code	CSB-YP856835DLU
Abbreviation	CG11412
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95RC0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MADAQAAAAG KKKYKNKKNS AEKNPNHNPN SSGQVEAQTP SNGHVQHQQE EATEDQEP AQ ELRGLLKKMH LCNGHGHKEQ EARPLGEVVN GHAHGHSSNN HIRCTSGSSN NNNSTHNNNS VDSSNNNRKQ RREGGDGGGS DSNSLKPEEK PITATSKTTA NIHPTTTTDP KPKVSEDAV EQGVHVATGS GHSREQERKQ PPSDYAEGAT PTITAQLQLP EPAISADEIV YKEYEAEHQM HDIMRLIQAE LSEPSIYTY RYFIYNWPKL CFLASHDNQY VGAIVCKLDM HMNVRRGYIA MLAVRKEYRK LKIGTTLVTK AIEAMLADNA DEVVLETEMR NQPALRLYEN LGFVRDKRLF RYYLNGVDAL RLKLWFR
Source	Yeast
Target Names	Naa30A
Protein Names	Recommended name: N-alpha-acetyltransferase 30 EC= 2.3.1.- Alternative name(s): N-acetyltransferase MAK3 homolog NatC catalytic subunit
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.