



Recombinant Human Protein FAM46B (FAM46B)

Product Code	CSB-MP856891HU
Abbreviation	FAM46B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96A09
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMPSESGAER RDRAAAQVGT AAATAVATAA PAGGGPDPEA LSAFPGRHLS GLSWPQVKRL DALLSEPIPI HGRGNFPTLS VQPRQIVQVV RSTLEEQGLH VHSVRLHGSA ASHVLHPESG LGYKDLDLVF RVDLRSEASF QLTKAVVLAC LLDFLPAGVS RAKITPLTLK EAYVQKLVKV CTDSDRWSLI SLSNKSGKNV ELKFVDSVRR QFEFSIDSFQ IILDSLFLFG QCSSTPMSEA FHPTVTGESL YGDFTEALEH LRHRVIATRS PEEIRGGGLL KYCHLLVRGF RPRPSTDVRA LQRYMCSRFF IDFPDLVEQR RTLERYLEAH FGGADAARRY ACLVTLHRVV NESTVCLMNH ERRQTLDLIA ALALQALAEQ GPAATAALAW RPPGTDGVVP ATVNYVTPV QPLLAHAYPT WLPCN
Source	Mammalian cell
Target Names	FAM46B
Protein Names	Recommended name: Protein FAM46B
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.