



Recombinant Human Pre-B-cell leukemia transcription factor-interacting protein 1 (PBXIP1)

Product Code	CSB-BP856896HU
Abbreviation	PBXIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96AQ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASCPDSDNS WVLAGESESLP VETLGPASRM DPESERALQA PHSPSKTDGK ELAGTMDGEG TLFQTESPQS GSILTEETEV KGTLEGDVCG VEPPGPGDTV VQGDLETTV VTGLPDTQD LEGQSPQSL PSTPKAAWIR EEGRCSSDD DTDVDMEGLR RRRGREAGPP QPMVPLAVEN QAGGEGAGGE LGISLNMCLL GALVLLGLGV LLFSGGLSES ETGPMEEVER QVLPDPEVLE AVGDRQDGLR EQLQAPVPPD SVPSLQNMGL LLDKLAKENQ DIRLLQAQLQ AQKEELQSLM HQPKGLEEN AQLRGALQQG EAFQRALESE LQQLRARLQG LEADCVRGPD GVCLSGGRGP QGDKAIREQG PREQEPELSF LKQKEQLEAE AQALRQELER QRLLGVSQQ DLERSLQDAS RGDPAHAGLA ELGHRLAQLK QGLENWGQDP GVSANASKAW HQKSHFQNSR EWSGKEKWWDD GQRDRKAEHW KHKKEESGRE RKKNWGGQED REPAGRWKEG RPRVEESGSK KEGKRQGPKE PPRKSGSFHS SGEKQKQPRW REGTKDSDHP LPSWAELLRP KYRAPQGCSG VDECARQEGL TFFGTELAPV RQQELASLLR TYLARLPWAG QLTKEPLSP AFFGEDGIFR HDRLRFRDFV DALEDSLEEV AVQQTGDDDE VDDFEDFIFS HFFGDKALKK RSGKKDKHSQ SPRAAGPREG HSHSHHHHR G
Source	Baculovirus
Target Names	PBXIP1
Protein Names	Recommended name: Pre-B-cell leukemia transcription factor-interacting protein 1 Alternative name(s): Hematopoietic PBX-interacting protein
Expression Region	1-731
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.