



Recombinant Human Dapper homolog 3 (DACT3)

| | |
|--------------------------|--|
| Product Code | CSB-YP856898HU |
| Abbreviation | DACT3 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q96B18 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MIRAFSFPVS PERGRLRGWL EGSLAGLCEL HWLRERQEYR VQQALRLAQP GMGGAEAEDE EDADEDEDAA AARRAAAAL EQLALPGLV WDLGQQLGDL SLESGGLEQE SGRSSGFYED PSSTGGPDSP PSTFCGDSGF SGSSSYGRLG PSEPRGIYAS ERPKSLGDAS PSAPEVVGAR AAVPRSFSAP YPTAGGSAGP EACSSAERRA RAGPFLTPSP LHAVAMRSPR PCGRPPTDSP DAGGAGRPLD GYISALLRRR RRRGAGQPR SPGGADGGPR RQNSVRQRPP DASPSPGSAR PAREPSLERV GGHPTSPAAL SRAWASSWES EAAPEPAAPP AAPSPDSPA EGRLVKAQYI PGAQAATRGL PGRAARRKPP PLTRGRSVEQ SPPRERPRAA GRRGRMAEAS GRRGSPRARK ASRSQSETSL LGRASAVPSG PPKYPTAERE EPRPPRPRRG PAPTAAQAA GSCRRWRSTA EIDAADGRRV RPRAPAARVP GPGPSPSAPQ RRLLYGCAGS DSECSAGRLG PLGRRGPAGG VGGGYGESES SASEGESPAF SSASSDSGDS GGLVWPQQLV AATAASGGGA GAGAPAGPAK VFKIKASHA LKKKILFRS GSLKVMTTV |
| Source | Yeast |
| Target Names | DACT3 |
| Protein Names | Recommended name: Dapper homolog 3 Alternative name(s): Antagonist of beta-catenin Dapper homolog 3 Arginine-rich region 1 protein |
| Expression Region | 1-629 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.