



Recombinant Human Protein FAM3D (FAM3D)

Product Code	CSB-YP856903HU
Abbreviation	FAM3D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BQ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	YMSFS MKTIRLPRWL AASPTKEIQV KKYKCGLIKP CPANYFAFKI CSGAANVVGP TMCFEDRMIM SPVKNNVGRG LNIALVNGTT GAVLGQKAFD MYSGDVMHLV KFLKEIPGGA LVLVASYDDP GTKMNDESRK LFSDLGSSYA KQLGFRDSWV FIGAKDLRGK SPFEQFLKNS PDTNKYEGWP ELLEMEGCMP PKPF
Source	Yeast
Target Names	FAM3D
Protein Names	Recommended name: Protein FAM3D
Expression Region	26-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.