



Recombinant Human Leukocyte receptor cluster member 1 (LENG1)

Product Code	CSB-EP856906HU-B
Abbreviation	LENG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BZ8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNILPKKSWH VRNKDNVARV RRDEAQAREE EKERERRVLL AQQEARTEFL RKKARHQNSL PELEAAEAGA PGSGPVDLFR ELLEEGKGV I RGNKEYEEEEK RQEKERQEKA LGILTYLGQS AAEAQTQPPW YQLPPGRGGP PPGPAPDEKI KSRLDPLREM QKHLGKKRQH GGDEGSRSRK EKEGSEKQRP KEPPSLDQLR AERLRREAAE RSRAEALLAR VQGRALQEGQ PEEDETDDRR RRYNSQFNPQ LARRPRQQDP HLTH
Source	E.coli
Target Names	LENG1
Protein Names	Recommended name: Leukocyte receptor cluster member 1
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.