



Recombinant Human ELL-associated factor 2 (EAF2)

Product Code	CSB-EP856910HU
Abbreviation	EAF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CJ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNSAAGFSLH DRRERVLKLG ESFEKQPRCA FHTVRYDFKP ASIDTSSEGY LEVGEQEQVT ITLPNIEGST PPVTVFKGSK KPYLKECILI INHDTGECRL EKLSSNITVK KTRVEGSSKI QYRKEQQQQQ MWNSARTPNL VKHSPSEDKM SPASPIDDIE RELKAEASLM DQMSSCDSSS DSKSSSSSSSS EDSSSDSEDE DCKSSTSDTG NCVSGHPTMT QYRIPDIDAS HNRFRDNSGL LMNTLRNDLQ LSEGSDDSD
Source	E.coli
Target Names	EAF2
Protein Names	Recommended name: ELL-associated factor 2 Alternative name(s): Testosterone-regulated apoptosis inducer and tumor suppressor protein
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.