



Recombinant Human FAS-associated factor 2 (FAF2)

Product Code	CSB-EP856912HU-B
Abbreviation	FAF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CS3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AAPEERDLT QEQTEKLLQF QDLTGIESMD QCRHTLEQHN WNIEAAVQDR LNEQEGVPSV FNPPSRPLQ VNTADHRIYS YVSRPQPRG LLGWGYLIM LPFRFTYYTI LDIFRFALRF IRPDPRSRVT DPVGDIVSFM HSFEKEYGRA HPVIFYQGTYS QALNDAKREL RFLLVYLHGD DHQDSDEFER NTLCAPEVIS LINTRMLFWA CSTNKPEGYR VSQALRENTY PFLAMIMLKD RRMTVVGRL GLIQPDDLIN QLTFIMDANQ TYLVSERLER EERNQTQVLR QQQDEAYLAS LRADQEKERK KREERERKRR KEEEVQQQKL AEERRRQNLQ EEKERKLECL PPEPSPDDPE SVKIIFKLPN DSRVERRFHF SQSLTVIHDF LFSLKESPEK FQIEANFPRR VLPCIPSEEW PNPPTLQEAG LSHTEVLFVQ DLTDE
Source	E.coli
Target Names	FAF2
Protein Names	Recommended name: FAS-associated factor 2 Alternative name(s): Protein ETEA UBX domain-containing protein 3B UBX domain-containing protein 8
Expression Region	2-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is highly expressed in peripheral blood of patients with atopic dermatitis (AD), compared to normal individuals. It may play a role in regulating the resistance to apoptosis that is observed in T cells and eosinophils of AD patients.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.